

Product datasheet for **MP201226**

Bhlhe41 Mouse qPCR Primer Pair (NM_024469)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79362
Forward Sequence:	GGCATTGGAGAAAGCAGTAGTC
Reverse Sequence:	AGTGGAACGCATCCAAGTCGGC
ACCN:	<u>NM_024469</u> , <u>NM_024469.1</u> , <u>NM_024469.2</u> , <u>BC023231</u> , <u>BC037742</u> , <u>BC051555</u> , <u>BC172625</u> , <u>BC172727</u> , <u>BE980915</u> , <u>BE985179</u>
UniProt ID:	<u>Q99PV5</u>
Synonyms:	6430520M22Rik; Bhlhb; Bhlhb2l; Bhlhb3; DEC2; Sh; Sharp1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)