

Product datasheet for **MP201246**

Bcl11a Mouse qPCR Primer Pair (NM_016707)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	14025
Forward Sequence:	CGGGATGAGTGCAGAATATGCC
Reverse Sequence:	GCAAGAGAAACCATGCACTGGTG
ACCN:	<u>BC010585</u> , <u>NM_016707</u> , <u>NM_016707.1</u> , <u>NM_016707.2</u> , <u>NM_016707.3</u> , <u>BB649527</u> , <u>BB662872</u> , <u>BB853515</u> , <u>BC051418</u> , <u>BF731393</u> , <u>BY735002</u>
UniProt ID:	<u>Q9QYE3</u>
Synonyms:	2810047E18Rik; BCL-11A; Ct1p1; D930021L15Rik; Evi9; Evi9a; Evi9b; Evi9c; mKIAA1809
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)