

Product datasheet for **MP201301**

Tspoap1 Mouse qPCR Primer Pair (NM_172449)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	207777
Forward Sequence:	GCAGTTCACCACACTTCCAGCA
Reverse Sequence:	ATGGTGACAGGGAGCCAAGTGA
ACCN:	BC066146 , NM_172449 , NM_172449.1 , NM_172449.2 , BC066146.1 , BC055691 , BC083181 , BQ769688 , BY233070
UniProt ID:	Q7TNF8
Synonyms:	Bzrap1; D230016K05; mKIAA0612; PRAX-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)