

## Product datasheet for **MP201332**

### Bfar Mouse qPCR Primer Pair (NM\_025976)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67118
Forward Sequence:	GAGTGAAGCCACCGCAGAATCT
Reverse Sequence:	GTATAGCAAGCCAAGTCTGGGTG
ACCN:	<u><a href="#">BC027221</a></u> , <u><a href="#">NM_025976</a></u> , <u><a href="#">NM_025976.1</a></u> , <u><a href="#">NM_025976.2</a></u> , <u><a href="#">NM_025976.3</a></u> , <u><a href="#">NM_025976.4</a></u> , <u><a href="#">NM_025976.5</a></u> , <u><a href="#">BB613414</a></u> , <u><a href="#">BY329099</a></u> , <u><a href="#">BY777666</a></u>
UniProt ID:	<u><a href="#">Q8R079</a></u>
Synonyms:	3010001A07Rik; 3110001I22Rik; AI666707; AW107665; Bar; Rnf47
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)