

Product datasheet for **MP201369**

Bmp15 Mouse qPCR Primer Pair (NM_009757)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12155
Forward Sequence:	GATTGGAGCGAAAATGGTGAGGC
Reverse Sequence:	GCTACCTGGTTTGATGCTAGAGG
ACCN:	<u>BC055363</u> , <u>NM_009757</u> , <u>NM_009757.1</u> , <u>NM_009757.2</u> , <u>NM_009757.3</u> , <u>NM_009757.4</u> , <u>NM_009757.5</u>
UniProt ID:	<u>Q9Z0L4</u>
Synonyms:	AU015375; AU018861; AU021453; Bmp-1; Bmp-15; C86824; C87336; GDF-9B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)