

Product datasheet for **MP201380**

Bmpr2 Mouse qPCR Primer Pair (NM_007561)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12168
Forward Sequence:	AGAGACCCAAGTCCCAGAAGC
Reverse Sequence:	TCTCCTCAGCACACTGTGCAGT
ACCN:	<u>NM_007561</u> , <u>NM_007561.1</u> , <u>NM_007561.2</u> , <u>NM_007561.3</u> , <u>NM_007561.4</u> , <u>BC148339</u> , <u>BC156868</u> , <u>BM940059</u>
UniProt ID:	<u>Q35607</u>
Synonyms:	2610024H22Rik; AL117858; AW546137; BB189135; BM; BMP; BMP-; BMP-2; BMPR-2; BMPR-II; BMPRII; BRK-3; Gm20272
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)