

Product datasheet for **MP201415**

Brip1 Mouse qPCR Primer Pair (NM_178309)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	237911
Forward Sequence:	CAGTAGCAAGCCACGGTGTTC
Reverse Sequence:	GCCAGTTGACTGGCTCTGTTG
ACCN:	<u>BC094252</u> , <u>NM_178309</u> , <u>NM_178309.1</u> , <u>NM_178309.2</u>
UniProt ID:	<u>Q5SXJ3</u>
Synonyms:	3110009N10Rik; 8030460J03Rik; BACH; Bach1; FACJ; Fancj; OF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)