

## Product datasheet for **MP201427**

### **Bcas3 Mouse qPCR Primer Pair (NM\_138681)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	192197
<b>Forward Sequence:</b>	GACTTCCAGGTCGTGGTTTGCA
<b>Reverse Sequence:</b>	ATGTGCTCCACTAAGGTCCCGT
<b>ACCN:</b>	<a href="#">BC118965</a> , <a href="#">BC118966</a> , <a href="#">NM_138681</a> , <a href="#">NM_138681.1</a> , <a href="#">NM_138681.2</a> , <a href="#">NM_138681.3</a> , <a href="#">NM_138681.4</a> , <a href="#">BC028339</a> , <a href="#">BC038173</a> , <a href="#">BC050233</a> , <a href="#">BC060690</a> , <a href="#">BC080775</a> , <a href="#">BC086693</a> , <a href="#">BY116917</a>
<b>UniProt ID:</b>	<a href="#">Q8CCN5</a>
<b>Synonyms:</b>	1500019F07Rik; 2610028P08Rik; AU021018; BC028339; K20D4; rudhira
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)