

Product datasheet for **MP201438**

Baz1b Mouse qPCR Primer Pair (NM_011714)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22385
Forward Sequence:	GTGCTGTCAAGAGCAGAAGGCT
Reverse Sequence:	CAGCCTTGTGTTCCGCATTCG
ACCN:	<u>NM_011714</u> , <u>NM_011714.1</u> , <u>NM_011714.2</u> , <u>BC141399</u> , <u>BC014741</u> , <u>BC171975</u> , <u>BF463363</u> , <u>BM114461</u>
UniProt ID:	<u>Q9Z277</u>
Synonyms:	C87820; Wbscr9; WSTF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)