

Product datasheet for **MP201441**

Brpf1 Mouse qPCR Primer Pair (NM_030178)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	78783
Forward Sequence:	CCAAGAGAAGGACACAGGCAAC
Reverse Sequence:	GTAGCGGTAAGCCTCCAAGTTC
ACCN:	BC046521 , NM_030178 , NM_030178.1 , NM_030178.2 , BC005647 , BC031594 , BC051640 , BC138361 , BC145256 , BY469507
UniProt ID:	B2RRD7
Synonyms:	4833438B11Rik; 4930540D11Rik; Brpf2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)