

Product datasheet for **MP201447**

Brd2 Mouse qPCR Primer Pair (NM_010238)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	14312
Forward Sequence:	CAATCCTCCAGACCACGATGTTG
Reverse Sequence:	TGGTAGAGGTCCTGGTTCCAGT
ACCN:	<u>NM_010238</u> , <u>NM_010238.1</u> , <u>NM_010238.2</u> , <u>NM_010238.3</u> , <u>BC138656</u> , <u>BC046809</u> , <u>BC085251</u> , <u>BC138655</u> , <u>BM947875</u> , <u>BY244178</u>
UniProt ID:	<u>Q7IJ13</u>
Synonyms:	AW228947; D17H6S113E; Frg-1; Fsrg-1; Fsrg1; mKIAA4005; Nat; Ring3; Rnf3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)