

## Product datasheet for **MP201474**

### Banp Mouse qPCR Primer Pair (NM\_016812)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	53325
Forward Sequence:	GGAAGACCCATTGAGCAACAGG
Reverse Sequence:	TTCGTTGTGGCTGTCGTCCTGA
ACCN:	<a href="#">BC021650</a> , <a href="#">BC022168</a> , <a href="#">BC062641</a> , <a href="#">NM_016812</a> , <a href="#">NM_016812.1</a> , <a href="#">NM_016812.2</a> , <a href="#">NM_016812.3</a> , <a href="#">NM_016812.4</a> , <a href="#">BB857942</a> , <a href="#">BC013339</a> , <a href="#">BU703599</a>
UniProt ID:	<a href="#">Q8VBU8</a>
Synonyms:	AA408158; SMAR1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)