

Product datasheet for **MP201525**

Cdh5 Mouse qPCR Primer Pair (NM_009868)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12562
Forward Sequence:	GAACGAGGACAGCAACTTCACC
Reverse Sequence:	GTTAGCGTGCTGGTTCCAGTCA
ACCN:	<u>BC054790</u> , <u>NM_009868</u> , <u>NM_009868.1</u> , <u>NM_009868.2</u> , <u>NM_009868.4</u> , <u>BP768525</u> , <u>BY078746</u>
UniProt ID:	<u>P55284</u>
Synonyms:	7B4; AA408225; Cd144; VE; VE-C; VE-Cad; Vec; VEcad; VECD
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)