

Product datasheet for **MP201543**

Cib1 Mouse qPCR Primer Pair (NM_011870)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23991
Forward Sequence:	TCACCTACCAGAGACAGCCTGA
Reverse Sequence:	TCACAAGCTGGCTCAGGTCTTC
ACCN:	<u>BC003714</u> , <u>BC054385</u> , <u>NM_011870</u> , <u>NM_011870.1</u> , <u>NM_011870.2</u> , <u>NM_011870.3</u> , <u>NM_011870.4</u> , <u>NM_011870.5</u> , <u>BB843350</u> , <u>BP770454</u>
UniProt ID:	<u>Q9Z0F4</u>
Synonyms:	Cibkip; Kip; Prkdcip
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)