

Product datasheet for **MP201562**

Cacnb2 Mouse qPCR Primer Pair (NM_023116)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12296
Forward Sequence:	CAACTCAGTAAGACCTCTTTGGC
Reverse Sequence:	CTGCTACCATCTGGACATTGAGG
ACCN:	BC109156 , BC109157 , NM_023116 , NM_023116.1 , NM_023116.2 , NM_023116.3 , NM_023116.4 , BC048944 , BC056183
UniProt ID:	Q8CC27
Synonyms:	AW060387; CAB2; Cavbeta2; Cchb2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)