

Product datasheet for **MP201580**

Cacna1h Mouse qPCR Primer Pair (NM_021415)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	58226
Forward Sequence:	GCACAAAGTGCTGGAGCCCTAT
Reverse Sequence:	GTGTGCGATGACTTTCTGGCAG
ACCN:	<u>NM_021415</u> , <u>NM_021415.2</u> , <u>NM_021415.3</u> , <u>NM_021415.4</u> , <u>BC138026</u> , <u>BC058206</u> , <u>BU612176</u> , <u>BU612210</u>
Synonyms:	alpha13.2; Cav3; Cav3.2; MNCb-1209
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)