

Product datasheet for **MP201585**

Camk1 Mouse qPCR Primer Pair (NM_133926)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	52163
Forward Sequence:	GGAGAACGAGATTGCCGTCTTAC
Reverse Sequence:	GGTCAAACAGCTCTCCACCTGA
ACCN:	<u>BC014825</u> , <u>NM_133926</u> , <u>NM_133926.1</u> , <u>NM_133926.2</u> , <u>BC025510</u> , <u>BC042494</u> , <u>BC057155</u> , <u>BC069882</u>
UniProt ID:	<u>Q91YS8</u>
Synonyms:	AI505105; CaMKIalpha; D6Ertd263e
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)