

Product datasheet for **MP201608**

Calm2 Mouse qPCR Primer Pair (NM_007589)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12314
Forward Sequence:	GATTAGAGAAGCATTCCGTGTGT
Reverse Sequence:	CCGTCAATGTCTGCTCCCTGA
ACCN:	BC021347 , BC051444 , BC100301 , NM_007589 , NM_007589.1 , NM_007589.2 , NM_007589.3 , NM_007589.4 , NM_007589.5 , BC010730 , BY095488 , BY660536
UniProt ID:	P0DP26
Synonyms:	1500001E21Rik; AL024017; Cam2; CamC
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)