

Product datasheet for **MP201670**

Cnrip1 Mouse qPCR Primer Pair (NM_029861)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	380686
Forward Sequence:	CGGCATCTATGACACAGAAGGTG
Reverse Sequence:	GGCAATGGTCTCGCTTGTGGTA
ACCN:	<u>BC087946, NM_029861, NM_029861.1, NM_029861.2, NM_029861.3</u>
UniProt ID:	<u>Q5M8N0</u>
Synonyms:	1500041B16Rik; 3110054C06Rik; 5330437A18Rik; AI854501; C2orf32
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)