

## Product datasheet for **MP201695**

### Chst14 Mouse qPCR Primer Pair (NM\_028117)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	72136
Forward Sequence:	CGCTTCCTCTACTGCTATGTCC
Reverse Sequence:	GTCTGCCAGAAACACCAAGTCAC
ACCN:	<u><a href="#">BC043700</a></u> , <u><a href="#">BC085479</a></u> , <u><a href="#">NM_028117</a></u> , <u><a href="#">NM_028117.1</a></u> , <u><a href="#">NM_028117.2</a></u> , <u><a href="#">NM_028117.3</a></u> , <u><a href="#">BC026886</a></u>
UniProt ID:	<u><a href="#">Q80V53</a></u>
Synonyms:	2600016L03Rik; D4ST-1; D4st1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)