

Product datasheet for **MP201775**

Csnk1g2 Mouse qPCR Primer Pair (NM_134002)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	103236
Forward Sequence:	CCAAACTCCGAGTCGGCAAGA
Reverse Sequence:	GTAAAAGCGGTACTCCAGGTGC
ACCN:	BC004839 , BC005750 , NM_134002 , NM_134002.1 , NM_134002.2 , BC035440 , BC063083 , BC128313
UniProt ID:	Q8BVP5
Synonyms:	2810429I12Rik; AI463719
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)