

Product datasheet for **MP201794**

Casp3 Mouse qPCR Primer Pair (NM_009810)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12367
Forward Sequence:	GGAGTCTGACTGGAAAGCCGAA
Reverse Sequence:	CTTCTGGCAAGCCATCTCCTCA
ACCN:	BC038825 , NM_009810 , NM_009810.1 , NM_009810.2 , NM_009810.3 , BC021500
UniProt ID:	P70677
Synonyms:	A830040C14Rik; AC-; AC-3; Casp; CASP-3; Caspase-3; CC3; CPP; CPP-32; CPP32; Lice; mld; mldy; SCA-1; Ya; Yama
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)