

## Product datasheet for **MP201807**

### Cat Mouse qPCR Primer Pair (NM\_009804)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12359
Forward Sequence:	CGGCACATGAATGGCTATGGATC
Reverse Sequence:	AAGCCTTCCTGCCTCTCCAACA
ACCN:	<u><a href="#">BC013447</a></u> , <u><a href="#">NM_009804</a></u> , <u><a href="#">NM_009804.1</a></u> , <u><a href="#">NM_009804.2</a></u> , <u><a href="#">BC047126</a></u> , <u><a href="#">BY093963</a></u>
UniProt ID:	<u><a href="#">P24270</a></u>
Synonyms:	2210418N07; Cas-1; Cas1; Cs-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)