

Product datasheet for **MP201815**

Ctnnd1 Mouse qPCR Primer Pair (NM_007615)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12388
Forward Sequence:	CTGTGATGGTGTTCCTGCTCTG
Reverse Sequence:	TGGGATGAGAGATTCCACAGGG
ACCN:	BC043108 , BC046589 , BC054544 , NM_007615 , NM_007615.1 , NM_007615.2 , NM_007615.3 , NM_007615.4 , BY061937 , BY422861
UniProt ID:	P30999
Synonyms:	AA409437; AU019353; Catns; Ctnnd; mKIAA0384; P120
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)