

Product datasheet for **MP201816**

Ctnnd2 Mouse qPCR Primer Pair (NM_008729)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18163
Forward Sequence:	GTGAGACTCCTTCGCAAGACCA
Reverse Sequence:	GGATAATCACCGCATTGGTCAGC
ACCN:	NM_008729 , NM_008729.1 , NM_008729.2 , BC138260 , BC138261 , BC145217 , BE847220 , BU611309 , NM_008729.3
UniProt ID:	O35927
Synonyms:	Catnd2; neurojugin; Nprap
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)