

Product datasheet for **MP201826**

Ctsc Mouse qPCR Primer Pair (NM_009982)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	13032
Forward Sequence:	TGCAGCACATCTTGAGGTCTC
Reverse Sequence:	GCAGTCCAAGACTTCTGAACGG
ACCN:	<u>BC067063</u> , <u>NM_009982</u> , <u>NM_009982.1</u> , <u>NM_009982.2</u> , <u>NM_009982.3</u> , <u>NM_009982.4</u> , <u>NM_009982.5</u> , <u>BQ561752</u> , <u>BY381435</u>
UniProt ID:	<u>P97821</u>
Synonyms:	AI047818; CatC; D; DP; DPP1; DPPI
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)