

Product datasheet for **MP201844**

Catsper3 Mouse qPCR Primer Pair (NM_029772)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	76856
Forward Sequence:	GGAGGTCAACAGGCTGATGAAC
Reverse Sequence:	GCTGAAGTCATCCAGAACCATGG
ACCN:	<u>BC089518</u> , <u>NM_029772</u> , <u>NM_029772.1</u> , <u>NM_029772.2</u> , <u>NM_029772.3</u> , <u>NM_029772.4</u>
UniProt ID:	<u>Q80W99</u>
Synonyms:	4921522D01Rik; CACRC
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)