

## Product datasheet for **MP201855**

### Cebpa Mouse qPCR Primer Pair (NM\_007678)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12606
Forward Sequence:	GCAAAGCCAAGAAGTCGGTGGA
Reverse Sequence:	CCTTCTGTTGCGTCTCCAGTT
ACCN:	<a href="#">BC011118</a> , <a href="#">BC028890</a> , <a href="#">BC051102</a> , <a href="#">BC058161</a> , <a href="#">NM_007678</a> , <a href="#">NM_007678.1</a> , <a href="#">NM_007678.2</a> , <a href="#">NM_007678.3</a> , <a href="#">BC106174</a>
UniProt ID:	<a href="#">P53566</a>
Synonyms:	C/ebp; C/ebpalpha; CBF-A; Ceb; Cebp
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)