

Product datasheet for **MP201859**

Cebpg Mouse qPCR Primer Pair (NM_009884)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12611
Forward Sequence:	GCTTACAGCAGGTTCTCAGCT
Reverse Sequence:	GGCGGTATTCGTCCTATTCCG
ACCN:	<u>BC011319</u> , <u>NM_009884</u> , <u>NM_009884.1</u> , <u>NM_009884.2</u> , <u>NM_009884.3</u> , <u>BE987706</u>
UniProt ID:	<u>P53568</u>
Synonyms:	C/EBPgamma; C/EBP[g]; GPE1-BP; Gpe1bp; Ig/EBP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

