

## Product datasheet for **MP201914**

### Cd300lf Mouse qPCR Primer Pair (NM\_145634)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	246746
Forward Sequence:	GATGCTGGCATTACTGGTGTGG
Reverse Sequence:	GGTTGTCAGTGAAGATGGTGG
ACCN:	<u><a href="#">BC057864</a>, <a href="#">NM_145634</a>, <a href="#">NM_145634.1</a>, <a href="#">NM_145634.2</a>, <a href="#">NM_145634.3</a>, <a href="#">BM227747</a></u>
UniProt ID:	<u><a href="#">Q6SIQ7</a></u>
Synonyms:	CLIM1; CLM-1; CLM1; Digr2; F730004D16Rik; IgSF13; IREM1; LMIR3; Pigr3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)