

## Product datasheet for **MP201922**

### Cd34 Mouse qPCR Primer Pair (NM\_133654)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12490
Forward Sequence:	AGGACAGCAGTAAGACCACACC
Reverse Sequence:	GTGTGGAGTTCCAGAGCCTGAA
ACCN:	<u><a href="#">BC006607</a></u> , <u><a href="#">NM_133654</a></u> , <u><a href="#">NM_133654.1</a></u> , <u><a href="#">NM_133654.2</a></u> , <u><a href="#">NM_133654.3</a></u> , <u><a href="#">BC066820</a></u> , <u><a href="#">BU058977</a></u> , <u><a href="#">BY249301</a></u> , <u><a href="#">BY327943</a></u>
UniProt ID:	<u><a href="#">Q64314</a></u>
Synonyms:	AU040960
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)