

Product datasheet for **MP201923**

Cd36 Mouse qPCR Primer Pair (NM_007643)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12491
Forward Sequence:	GGACATTGAGATTCTTTTCCTCTG
Reverse Sequence:	GCAAAGGCATTGGCTGGAAGAAC
ACCN:	<u>BC010262</u> , <u>NM_007643</u> , <u>NM_007643.1</u> , <u>NM_007643.2</u> , <u>NM_007643.3</u> , <u>NM_007643.4</u> , <u>BB534670</u> , <u>BC027691</u> , <u>BQ926044</u> , <u>BY126014</u>
UniProt ID:	<u>Q08857</u>
Synonyms:	FAT; GPIV; Scarb3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)