

Product datasheet for **MP201929**

Cd40lg Mouse qPCR Primer Pair (NM_011616)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	21947
Forward Sequence:	GAACTGTGAGGAGATGAGAAGGC
Reverse Sequence:	TGGCTTCGCTTACAACGTGTGC
ACCN:	BC119225 , NM_011616 , NM_011616.1 , NM_011616.2
UniProt ID:	P27548
Synonyms:	CD40-L; Cd40l; CD154; gp39; HIGM1; IGM; IMD3; Ly-62; Ly62; T-BAM; Tnfsf5; Tnlg8b; TRAP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)