

Product datasheet for **MP201956**

Cd86 Mouse qPCR Primer Pair (NM_019388)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12524
Forward Sequence:	ACGTATTGGAAGGAGATTACAGCT
Reverse Sequence:	TCTGTCAGCGTTACTATCCCGC
ACCN:	<u>BC013807</u> , <u>NM_019388</u> , <u>NM_019388.1</u> , <u>NM_019388.2</u> , <u>NM_019388.3</u>
UniProt ID:	<u>P42082</u>
Synonyms:	B7; B7-2; B7.2; B70; Cd28l2; CLS1; ETC-1; Ly-58; Ly58; MB7; MB7-2; TS/A-2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)