

Product datasheet for **MP201965**

Cdc14b Mouse qPCR Primer Pair (NM_172587)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	218294
Forward Sequence:	GACACAAGGAGACAGACTTCGG
Reverse Sequence:	CACAATCCACCAGATGGCAACAC
ACCN:	<u>BC049794</u> , <u>BC057357</u> , <u>NM_172587</u> , <u>NM_172587.1</u> , <u>NM_172587.2</u> , <u>NM_172587.3</u> , <u>BY199064</u>
UniProt ID:	<u>Q6PFY9</u>
Synonyms:	2810432N10Rik; A530086E13Rik; AA472821; Cdc14B1; CDC14B3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)