

Product datasheet for **MP201985**

Cisd2 Mouse qPCR Primer Pair (NM_025902)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67006
Forward Sequence:	GATCTGTGTCTCACCAAAGCAGC
Reverse Sequence:	GGATGAGAGGACCCACGTTATC
ACCN:	BC058279 , NM_025902 , NM_025902.1 , NM_025902.2 , NM_025902.3 , BC025614 , BC037182 , BC052543 , BE944243
UniProt ID:	Q9CQB5
Synonyms:	1500009M05Rik; 1500026J14Rik; 1500031D15Rik; AI848398; B630006A20Rik; Miner1; Noxp70
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)