

Product datasheet for **MP202029**

Cwc22 Mouse qPCR Primer Pair (NM_030560)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80744
Forward Sequence:	CCAGAGACAACCCACGCAATAC
Reverse Sequence:	TGGTTTCTGGGCTACGATGACC
ACCN:	<u>BC003993</u> , <u>NM_030560</u> , <u>NM_030560.1</u> , <u>NM_030560.2</u> , <u>NM_030560.3</u> , <u>NM_030560.4</u> , <u>NM_030560.5</u> , <u>BC003993.1</u> , <u>BC003273</u> , <u>BU613482</u>
UniProt ID:	<u>Q8C5N3</u>
Synonyms:	AA684037; AII73004; AL022752; B230213M24; mKIAA1604
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

