

Product datasheet for **MP202228**

Ceacam1 Mouse qPCR Primer Pair (NM_001039185)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26365
Forward Sequence:	CCTCTATTCCAGGAAGTCTGGC
Reverse Sequence:	GTTCAGGACAGTGTATGCGACG
ACCN:	<u>BC016891</u> , <u>NM_001039185</u> , <u>NM_001039185.1</u> , <u>BC114998</u> , <u>BC121189</u>
UniProt ID:	<u>P31809</u>
Synonyms:	bb-1; Bgp; Bgp1; C-CAM; Cc1; CD66a; Cea-1; Cea-7; Cea1; Cea7; Hv-2; Hv2; mCEA1; Mhv-1; MHVR; MHVR1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)