

Product datasheet for **MP202242**

Cadm3 Mouse qPCR Primer Pair (NM_053199)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	94332
Forward Sequence:	GTGCTCCATCTTCACTATGCCTG
Reverse Sequence:	AGAGTGGCTGTCTCCTTTTCCC
ACCN:	BC115417 , BC115418 , NM_053199 , NM_053199.1 , NM_053199.2 , NM_053199.3 , BC029659 , BC079837 , BY289385
UniProt ID:	Q99N28
Synonyms:	BlgR; Igsf4b; Necl-1; Necl1; SynCAM3; Tsl11
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)