

Product datasheet for **MP202260**

Cdc27 Mouse qPCR Primer Pair (NM_145436)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	217232
Forward Sequence:	AGTCGAGAGGTCACACCAGTAC
Reverse Sequence:	GTGAAAAGCGGGAAGTTCTCC
ACCN:	NM_145436 , NM_145436.1 , NM_145436.2 , BC023187 , BC023730 , BC040245 , BC157955 , BC160241 , BC172100 , BC172149 , BY013739 , BY036319
UniProt ID:	A2A6Q5
Synonyms:	AI452358; APC3; BC023187
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)