

Product datasheet for **MP202308**

Cenpe Mouse qPCR Primer Pair (NM_173762)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	229841
Forward Sequence:	AGGATCATGCCACCGAGAAGAC
Reverse Sequence:	GCTGTGTCTCTTGAGTTTCTGG
ACCN:	NM_173762 , NM_173762.1 , NM_173762.2 , NM_173762.3 , NM_173762.4 , BC039755 , BC049989 , BC052843 , BC059032 , BC080703 , BC106096
UniProt ID:	Q6RT24
Synonyms:	312kDa; AU019344; BC049989; C530022J18; CENP-E; Kif10
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)