

Product datasheet for **MP202310**

Cenph Mouse qPCR Primer Pair (NM_021886)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26886
Forward Sequence:	GGATGTTGACAAGATGGAGAATTC
Reverse Sequence:	CTGCCCAGTTGGTTTACTCGC
ACCN:	<u>BC025084, NM_021886, NM_021886.1</u>
UniProt ID:	<u>Q9QYM8</u>
Synonyms:	1700021111Rik; 2410018A12Rik; 2610042E16Rik; 2810046K12Rik; AU044255; CENP-H; ENP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

