

Product datasheet for **MP202327**

Cep57 Mouse qPCR Primer Pair (NM_026665)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	74360
Forward Sequence:	AACCCATGTTTCAGAGCCAGCTG
Reverse Sequence:	CCTTTCCTGCTCTTCTTCACGG
ACCN:	<u>BC050785</u> , <u>NM_026665</u> , <u>NM_026665.1</u> , <u>NM_026665.2</u> , <u>NM_026665.3</u> , <u>NM_026665.4</u> , <u>BC050785.1</u> , <u>BB846596</u> , <u>BC025588</u> , <u>BC030353</u> , <u>BC056968</u> , <u>BY134402</u>
UniProt ID:	<u>Q8CEE0</u>
Synonyms:	3110002L15Rik; 4921510P06Rik; 4931428M20Rik; AI467480; mKIAA0092; Tsp57
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

