

Product datasheet for **MP202336**

Cerk Mouse qPCR Primer Pair (NM_145475)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	223753
Forward Sequence:	GGCACAAACGACGCAGAGACAT
Reverse Sequence:	GAAACCGTAGCCCAGCAGAGAA
ACCN:	BC094253 , NM_145475 , NM_145475.1 , NM_145475.2 , NM_145475.3 , NM_145475.4 , BC038149 , BU756835
UniProt ID:	Q8K4Q7
Synonyms:	AI848610; D330016D08Rik; mKIAA1646
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)