

Product datasheet for **MP202370**

Xcl1 Mouse qPCR Primer Pair (NM_008510)

Product data:

| | |
|--------------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 16963 |
| Forward Sequence: | CAAGACCTATATCATCTGGGAGG |
| Reverse Sequence: | TGCCATCCACAGTCTTGATCGC |
| ACCN: | BC062249 , NM_008510 , NM_008510.1 , NM_008510.2 |
| UniProt ID: | P47993 |
| Synonyms: | AI661682; ATAC; Lptn; LTN; SCM-1; SCM-1a; Scyc1 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)