

Product datasheet for **MP202376**

Ccl2 Mouse qPCR Primer Pair (NM_011333)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20296
Forward Sequence:	GCTACAAGAGGATCACCAGCAG
Reverse Sequence:	GTCTGGACCCATTCTTCTGG
ACCN:	<u>NM_011333</u> , <u>NM_011333.1</u> , <u>NM_011333.2</u> , <u>NM_011333.3</u> , <u>BC145867</u> , <u>BC055070</u> , <u>BC145869</u>
UniProt ID:	<u>P10148</u>
Synonyms:	AI323594; HC11; JE; MCA; MCAF; MCP; MCP-; MCP-1; MCP1; Scy; Scya2; Sig; Sigje; SMC-C; SMC-CF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)