

## Product datasheet for **MP202379**

### Ccl21a Mouse qPCR Primer Pair (NM\_011124)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18829
Forward Sequence:	GGGTCAGGACTGCTGCCTTAAG
Reverse Sequence:	AGCTCAGGCTTAGAGTGCTTCC
ACCN:	<a href="#">BC025974</a> , <a href="#">BC028747</a> , <a href="#">BC038120</a> , <a href="#">NM_011124</a> , <a href="#">NM_011124.1</a> , <a href="#">NM_011124.2</a> , <a href="#">NM_011124.3</a> , <a href="#">NM_011124.4</a>
UniProt ID:	<a href="#">P84444</a>
Synonyms:	6CKBAC2; 6Ckine; ALP; AW987545; CKb9; plt; Scya21; SCYA21a; Scya21b; SLC; Tca4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)