

Product datasheet for **MP202383**

Ccl25 Mouse qPCR Primer Pair (NM_009138)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20300
Forward Sequence:	AAGGCTAGTCCACTGGAAGAGC
Reverse Sequence:	GTGGCACTCCTCACGCTTGAC
ACCN:	<u>BC117033</u> , <u>BC117037</u> , <u>NM_009138</u> , <u>NM_009138.1</u> , <u>NM_009138.2</u> , <u>NM_009138.3</u>
UniProt ID:	<u>O35903</u>
Synonyms:	A130072A22Rik; AI852536; CKb15; Scya25; TECK
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)