

Product datasheet for **MP202399**

Ccr4 Mouse qPCR Primer Pair (NM_009916)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12773
Forward Sequence:	GGACTAGGTCTGTGCAAGATCG
Reverse Sequence:	TGCCTTCAAGGAGAATACCGCG
ACCN:	<u>BC117041</u> , <u>BC119171</u> , <u>NM_009916</u> , <u>NM_009916.1</u> , <u>NM_009916.2</u>
UniProt ID:	<u>P51680</u>
Synonyms:	C-C CKR-4; CHEMR1; Cmkbr4; LESTR; Sdf1r
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)