

Product datasheet for **MP202409**

Cxcl1 Mouse qPCR Primer Pair (NM_008176)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	14825
Forward Sequence:	TCCAGAGCTTGAAGGTGTTGCC
Reverse Sequence:	AACCAAGGGAGCTTCAGGGTCA
ACCN:	<u>BC132502</u> , <u>BC132504</u> , <u>NM_008176</u> , <u>NM_008176.1</u> , <u>NM_008176.2</u> , <u>NM_008176.3</u> , <u>BC037997</u> , <u>BQ031102</u> , <u>BY335927</u>
UniProt ID:	<u>P12850</u>
Synonyms:	F; Fsp; Gr; gro; Gro1; KC; Mg; Mgsa; N5; N51; Scyb; Scyb1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)